

## Philips Respironics Trilogy 202

**Battery life approx. 3 hours**



**Two types of single limb circuit work with the Trilogy:**

Intersurgical **5804000**



Respironics **1065836**

Discard long thin tube  
and close port on circuit



**Expiratory Valve – DO NOT OCCLUDE**

1. Connect the circuit at side of the ventilator (see photo above). Connect O2 hose to port. Plug in.

2. Press and hold the ▼ and alarm silence buttons simultaneously for at least 5 seconds to enter **Setup** menu.

3. Use ▼ to scroll to **Settings and Alarms**. Press button under **Select**.

4. Set desired mode (**PC-SIMV, ST, CPAP**) and appropriate settings.

▲ ▼ to scroll to parameter or adjust, button under **Modify** to confirm to change, and button under **Ok** to confirm each setting.

Ensure **AUTO-TRAK** trigger (monitors and compensates for leaks) is **ON**. Note Circuit type = **passive** and dual prescription = **OFF**.

5. Set alarms (included in Settings and Alarms menu). **ENSURE** alarms for **Low MV** and **Circuit Disconnect** are **ON**.

6. Press button under **Finish**. If **mode changed** check details correct on screen, press button under **Yes** button then **Exit**. Mode not changed just press **Exit**. Ventilator turns off after **Exit** selected.

7. Turn on using **ON/OFF** button. Device performs a self-test and **VENTILATION WILL COMMENCE** with the previous settings.

**NOTE** device connects directly to oxygen outlet. Value displayed is set O<sub>2</sub>. Can deliver 100% O<sub>2</sub> for 2 mins via button under **100% O<sub>2</sub>** on monitoring page.

### Troubleshooting:

Circuit type **must be passive** – Setup Menu (step 2) -> Settings and alarms -> Circuit type **Passive**

Keypad locked – press and hold the right button for 5s to unlock. Turn function on/off in **Options**

Unable to view full details on display – change to detailed view under **Options**

If limited menu access – press ▼ and mute button together to re-access the menu and change to full under **Options**.

